




21 CFR Part 11 Compliance Readiness


User administration
(local or active directory authentication)



Comprehensive activity logging (by object, by user)

Date	Person	Location	Change	Comment
30.06.2021 16:02:37	Tester	Test	Imported from 'C:\Users\j\Downloads\Your Sample_01.spc' with LabCognition.IO.GalacticSpd&Import	
30.06.2021 16:02:37	Tester	Test	Label 'Filename' inserted with value 'Your Sample_01'	
30.06.2021 16:02:37	Tester	Test	Label 'DataSourceUrl' inserted with value 'C:\Users\j\Downloads\Your Sample_01.spc'	
30.06.2021 16:02:48	Tester	Test	Derivative: Order=1, Algorithm=Savitzky Golay, Smoothing Window=Nine Points	

Electronic signature & audit trail for data objects



Data and functional access permission schemes

Release and unrelease objects in keeping with organizational data control strategies